



# CSRNP3 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA855508NC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8WYN3
<b>Immunogen</b>	Peptide sequence from Human Cysteine/serine-rich nuclear protein 3 protein (518-536AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cysteine/serine-rich nuclear protein 3 (CSRNP-3) (Protein FAM130A2) (TGF-beta-induced apoptosis protein 2) (TAIP-2), CSRNP3, FAM130A2 TAIP2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	CSRNP3