



# CPNE1 Antibody, HRP conjugated

Product Code	CSB-PA859524LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99829
Immunogen	Recombinant Human Copine-1 protein (326-537AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Copine-1 (Chromobindin 17) (Copine I), CPNE1, CPN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CPNE1