



OSBPL7 Antibody, HRP conjugated

Product Code	CSB-PA861186LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZF2
Immunogen	Recombinant Human Oxysterol-binding protein-related protein 7 protein (1-249AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	autophagosome, cytosol, perinuclear endoplasmic reticulum, plasma membrane, cholesterol binding, sterol transporter activity, bile acid biosynthetic process, cellular response to cholesterol, positive regulation of proteasomal protein catabolic process, regulation of autophagy
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxysterol-binding protein-related protein 7 (ORP-7) (OSBP-related protein 7), OSBPL7, ORP7
Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OSBPL7