



OSBPL7 Antibody, Biotin conjugated

Product Code	CSB-PA861186LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZF2
Immunogen	Recombinant Human Oxysterol-binding protein-related protein 7 protein (1-249AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxysterol-binding protein-related protein 7 (ORP-7) (OSBP-related protein 7), OSBPL7, ORP7
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OSBPL7