



SNRK Antibody, FITC conjugated

Product Code	CSB-PA868301NC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRH2
Immunogen	Peptide sequence from Human SNF-related serine/threonine-protein kinase protein (340-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	May play a role in hematopoietic cell proliferation or differentiation. Potential mediator of neuronal apoptosis.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SNF-related serine/threonine-protein kinase (EC 2.7.11.1) (SNF1-related kinase), SNRK, KIAA0096 SNFRK
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SNRK