



CELSR1 Antibody, HRP conjugated

Product Code	CSB-PA868360NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYQ6
Immunogen	Peptide sequence from Human Cadherin EGF LAG seven-pass G-type receptor 1 protein (345-363AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin EGF LAG seven-pass G-type receptor 1 (Cadherin family member 9) (Flamingo homolog 2) (hFmi2), CELSR1, CDHF9 FMI2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CELSR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.