



CELSR1 Antibody, Biotin conjugated

Product Code	CSB-PA868360ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYQ6
Immunogen	Peptide sequence from Human Cadherin EGF LAG seven-pass G-type receptor 1 protein (345-363AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin EGF LAG seven-pass G-type receptor 1 (Cadherin family member 9) (Flamingo homolog 2) (hFmi2), CELSR1, CDHF9 FMI2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CELSR1