



B4GALT6 Antibody, Biotin conjugated

Product Code	CSB-PA871391LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBX8
Immunogen	Recombinant Human Beta-1,4-galactosyltransferase 6 protein (70-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-1,4-galactosyltransferase 6 (Beta-1,4-GalTase 6) (Beta4Gal-T6) (b4Gal-T6) (UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 6) (UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 6) [Includes: Glucosylceramide beta-1,4-galactosyltransferase (EC 2.4.1.274) (Lactosylceramide synthase) (LacCer synthase) (UDP-Gal:glucosylceramide beta-1,4-galactosyltransferase)], B4GALT6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	B4GALT6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.