



DUSP13 Antibody, Biotin conjugated

Product Code	CSB-PA883436LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UII6
Immunogen	Recombinant Human Dual specificity protein phosphatase 13 isoform B protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 13 isoform B (DUSP13B) (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity phosphatase SKRP4) (Testis- and skeletal-muscle-specific DSP), DUSP13, DUSP13B TMDP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DUSP13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.