



TBK1 Antibody, HRP conjugated

Product Code	CSB-PA890690NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHD2
Immunogen	Peptide sequence from Human Serine/threonine-protein kinase TBK1 protein (644-661AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase TBK1 (EC 2.7.11.1) (NF-kappa-B-activating kinase) (T2K) (TANK-binding kinase 1), TBK1, NAK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.