



FEM1B Antibody, FITC conjugated

Product Code	CSB-PA891718LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK73
Immunogen	Recombinant Human Protein fem-1 homolog B protein (400-627AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Component of an E3 ubiquitin-protein ligase complex, in which it may act as a substrate recognition subunit. Involved in apoptosis by acting as a death receptor-associated protein that mediates apoptosis. Also involved in glucose homeostasis in pancreatic islet. Functions as an adapter/mediator in replication stress-induced signaling that leads to the activation of CHEK1.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein fem-1 homolog B (FEM1b) (FEM1-beta) (Fem-1-like death receptor-binding protein alpha) (Fem-1-like in apoptotic pathway protein alpha) (F1A-alpha), FEM1B, F1AA KIAA0396
Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FEM1B