



SLC39A10 Antibody, Biotin conjugated

Product Code	CSB-PA891736LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULF5
Immunogen	Recombinant Human Zinc transporter ZIP10 protein (41-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter ZIP10 (Solute carrier family 39 member 10) (Zrt- and Irt-like protein 10) (ZIP-10), SLC39A10, KIAA1265 ZIP10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC39A10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.