



TENM1 Antibody, Biotin conjugated

Product Code	CSB-PA892339LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKZ4
Immunogen	Recombinant Human Teneurin-1 protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Teneurin-1 (Ten-1) (Protein Odd Oz/ten-m homolog 1) (Tenascin-M1) (Ten-m1) (Teneurin transmembrane protein 1) [Cleaved into: Ten-1 intracellular domain (IDten-1) (Ten-1 ICD); Teneurin C-terminal-associated peptide (TCPA-1) (Ten-1 extracellular domain) (Ten-1 ECD)], TENM1, ODZ1 TNM1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TENM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.