



NENF Antibody, Biotin conjugated

Product Code	CSB-PA892354ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMX5
Immunogen	Peptide sequence from Human Neudesin protein (149-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neudesin (Cell immortalization-related protein 2) (Neuron-derived neurotrophic factor) (Protein GIG47) (Secreted protein of unknown function) (SPUF protein), NENF, CIR2 SPUF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NENF