



NENF Antibody, Biotin conjugated

Product Code	CSB-PA892354ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMX5
Immunogen	Peptide sequence from Human Neudesin protein (149-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neudesin (Cell immortalization-related protein 2) (Neuron-derived neurotrophic factor) (Protein GIG47) (Secreted protein of unknown function) (SPUF protein), NENF, CIR2 SPUF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NENF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.