



IVNS1ABP Antibody, HRP conjugated

Product Code	CSB-PA896937LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6Y0
Immunogen	Recombinant Human Influenza virus NS1A-binding protein (401-642AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Influenza virus NS1A-binding protein (NS1-BP) (NS1-binding protein) (Aryl hydrocarbon receptor-associated protein 3), IVNS1ABP, ARA3 FLARA3 KIAA0850 NS1 NS1BP
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	IVNS1ABP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.