



# IVNS1ABP Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA896937LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y6Y0
<b>Immunogen</b>	Recombinant Human Influenza virus NS1A-binding protein (401-642AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Influenza virus NS1A-binding protein (NS1-BP) (NS1-binding protein) (Aryl hydrocarbon receptor-associated protein 3), IVNS1ABP, ARA3 FLARA3 KIAA0850 NS1 NS1BP
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Microbiology
<b>Target Names</b>	IVNS1ABP
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.