



ARL2BP Antibody, HRP conjugated

Product Code	CSB-PA897097HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2Y0
Immunogen	Recombinant Human ADP-ribosylation factor-like protein 2-binding protein (1-163AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-like protein 2-binding protein (ARF-like 2-binding protein) (ARL2-binding protein) (Binder of ARF2 protein 1), ARL2BP, BART BART1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARL2BP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.