



ARL2BP Antibody, HRP conjugated

Product Code	CSB-PA897097HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2Y0
Immunogen	Recombinant Human ADP-ribosylation factor-like protein 2-binding protein (1-163AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Together with ARL2, plays a role in the nuclear translocation, retention and transcriptional activity of STAT3. May play a role as an effector of ARL2.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-like protein 2-binding protein (ARF-like 2-binding protein) (ARL2-binding protein) (Binder of ARF2 protein 1), ARL2BP, BART BART1
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARL2BP