




FEN1 Antibody

Product Code	CSB-PA008585LA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O65251
Immunogen	Recombinant Arabidopsis thaliana Flap endonuclease 1 protein (85-383AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flap endonuclease 1, FEN-1, Flap structure-specific endonuclease 1, FEN1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FEN1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: FEN1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 37 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.