



Slfn9 Antibody, HRP conjugated

Product Code	CSB-PA106533LB21MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1ARD6
Immunogen	Recombinant Mouse Schlafen family member 9 protein (1-83AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Schlafen family member 9, Schlafen-9, Slfn9
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Slfn9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.