



Met Antibody, FITC conjugated

Product Code	CSB-PA27899C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	V5KYY7
Immunogen	Recombinant Locusta migratoria Met protein (1-329AA)
Raised In	Rabbit
Species Reactivity	Locusta migratoria
Tested Applications	ELISA
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Met protein, Met
Species	Locusta migratoria (Migratory locust)
Research Area	Signal Transduction
Target Names	Met