



Met Antibody, Biotin conjugated

Product Code	CSB-PA27899D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	V5KYY7
Immunogen	Recombinant Locusta migratoria Met protein (1-329AA)
Raised In	Rabbit
Species Reactivity	Locusta migratoria
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Met protein, Met
Immunogen Species	Locusta migratoria (Migratory locust)
Research Area	Signal Transduction
Target Names	Met
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.