



crfb12 Antibody, Biotin conjugated

Product Code	CSB-PA29691D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1RM05
Immunogen	Peptide sequence from Zebrafish Cytokine receptor family member B12 protein (106-130AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytokine receptor family member B12, crfb12
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	crfb12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.