



ifnlr1 Antibody, HRP conjugated

Product Code	CSB-PA29701B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D9IL84
Immunogen	Peptide sequence from Zebrafish Cytokine receptor family member b14 protein (106-130AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytokine receptor family member b14, ifnlr1, crfb14
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	ifnlr1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.