



# ifnlr1 Antibody, HRP conjugated

Product Code	CSB-PA29701B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D9IL84
Immunogen	Peptide sequence from Zebrafish Cytokine receptor family member b14 protein (106-130AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytokine receptor family member b14, ifnlr1, crfb14
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	ifnlr1