



CD4L-2 Antibody

Product Code	CSB-PA29721A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A1B2FIB0
Immunogen	Peptide sequence from Grass carp CD4-like protein 2 protein (59-77AA)
Raised In	Rabbit
Species Reactivity	Grass carp
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD4L-2
Immunogen Species	Ctenopharyngodon idella (Grass carp) (Leuciscus idella)
Research Area	Others
Target Names	CD4L-2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.