



# CD4L-2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA29721C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	A0A1B2FIB0
<b>Immunogen</b>	Peptide sequence from Grass carp CD4-like protein 2 protein (59-77AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Grass carp
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	CD4L-2
<b>Immunogen Species</b>	Ctenopharyngodon idella (Grass carp) (Leuciscus idella)
<b>Research Area</b>	Others
<b>Target Names</b>	CD4L-2