



plac8.1 Antibody, Biotin conjugated

Product Code	CSB-PA29731D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q498P1
Immunogen	Peptide sequence from Zebrafish Placenta-specific 8, tandem duplicate 1 protein (66-82AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Placenta-specific 8, tandem duplicate 1, Zgc:114201, plac8.1, zgc:114201
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	plac8.1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.