



plac8.2 Antibody, FITC conjugated

Product Code	CSB-PA29741C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2VD52
Immunogen	Peptide sequence from Zebrafish Zgc:158845 protein (5-19AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zgc:158845 protein, plac8.2
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	plac8.2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.