



CUP1-1 Antibody, HRP conjugated

Product Code	CSB-PA317846DB01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CX80
Immunogen	Recombinant Saccharomyces cerevisiae Copper metallothionein 1-1 protein (9-61AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Copper metallothionein 1-1, Cu-MT, Cu-metallothionein, Copper chelatin, Copper thionein, CUP1-1, MTH1
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Others
Target Names	CUP1-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.