



pdeH Antibody, HRP conjugated

Product Code	CSB-PA336327LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37646
Immunogen	Recombinant Escherichia coli Cyclic di-GMP phosphodiesterase PdeH protein (1-255AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclic di-GMP phosphodiesterase PdeH, pdeH, yhjH
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	pdeH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.