




Copper-metallothionein Antibody

Product Code	CSB-PA345921DA01HDB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55947
Immunogen	Recombinant Roman snail Copper-metallothionein protein (1-64AA)
Raised In	Rabbit
Species Reactivity	Roman snail
Tested Applications	ELISA, WB; Recommended dilution: WB:1:2000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Copper-metallothionein, Cu-MT
Immunogen Species	Helix pomatia (Roman snail) (Edible snail)
Research Area	Others
Target Names	Copper-metallothionein

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Copper-metallothionein antibody at 1:2500

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 32.7 kDa
Observed band size: 35 kDa