



grxD Antibody, FITC conjugated

Product Code	CSB-PA359726LC01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AC69
Immunogen	Recombinant Escherichia coli Glutaredoxin 4 protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaredoxin 4, Grx4, Monothiol glutaredoxin, grxD, ydhD
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	grxD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.