




hns Antibody

Product Code	CSB-PA359764LA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0ACF8
Immunogen	Recombinant Escherichia coli DNA-binding protein H-NS protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein H-NS, Heat-stable nucleoid-structuring protein, Histone-like protein HLP-II, Protein B1, Protein H1, hns, bgIY, cur, drdX, hnsA, msyA, osmZ, pilG, topS
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	hns

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: hns antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 19 kDa

Observed band size: 21 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.