



hns Antibody, Biotin conjugated

Product Code	CSB-PA359764LD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0ACF8
Immunogen	Recombinant Escherichia coli DNA-binding protein H-NS protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein H-NS, Heat-stable nucleoid-structuring protein, Histone-like protein HLP-II, Protein B1, Protein H1, hns, bglY, cur, drdX, hnsA, msyA, osmZ, pilG, topS
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	hns