




nfuA Antibody

Product Code	CSB-PA540344LA01ENT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1X760
Immunogen	Recombinant Escherichia coli Fe/S biogenesis protein NfuA protein (1-191AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fe/S biogenesis protein NfuA, nfuA
Immunogen Species	Escherichia coli (strain K12 / DH10B)
Research Area	Others
Target Names	nfuA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: nfuA antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 25 kDa

Observed band size: 32 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.