



ACS8 Antibody, HRP conjugated

Product Code	CSB-PA871226LB01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9T065
Immunogen	Recombinant Arabidopsis thaliana 1-aminocyclopropane-1-carboxylate synthase 8 protein (1-469AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	1-aminocyclopropane-1-carboxylate synthase 8, ACC synthase 8, S-adenosyl-L-methionine methylthioadenosine-lyase 8, ACS8
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	ACS8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.