




TTYH2 Antibody

Product Code	CSB-PA871565LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSA4
Immunogen	Recombinant Human Protein tweety homolog 2 protein (455-534AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein tweety homolog 2, hTTY2, TTYH2, C17orf29
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTYH2

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: TTYH2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 22.9 kDa

Observed band size: 28 kDa