



TTYH2 Antibody, Biotin conjugated

Product Code	CSB-PA871565LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSA4
Immunogen	Recombinant Human Protein tweety homolog 2 protein (455-534AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein tweety homolog 2, hTTY2, TTYH2, C17orf29
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTYH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.