




grxD Antibody

Product Code	CSB-PA359726LA11ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AC69
Immunogen	Recombinant Escherichia coli Glutaredoxin 4 protein (24-115AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaredoxin 4, Grx4, Monothiol glutaredoxin, grxD, ydhD
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	grxD

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: grxD antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 37 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.