



nfuA Antibody, FITC conjugated

Product Code	CSB-PA540344EC11ENT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1X760
Immunogen	Recombinant Escherichia coli Fe/S biogenesis protein NfuA protein (31-128AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fe/S biogenesis protein NfuA, nfuA
Immunogen Species	Escherichia coli (strain K12 / DH10B)
Research Area	Others
Target Names	nfuA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.