



Allergen Ara h 1, clone P17 Antibody, HRP conjugated

Product Code	CSB-PA331779LB01ANE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43237
Immunogen	Recombinant Peanut Allergen Ara h 1, clone P17 protein (26-216AA)
Raised In	Rabbit
Species Reactivity	Peanut
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Allergen Ara h 1, clone P17 (Allergen Ara h I) (allergen Ara h 1)
Immunogen Species	Peanut
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.