




Recombinant Macaca fascicularis Tumor necrosis factor ligand superfamily member 6 (FASLG), partial

Product Code	CSB-EP008434MOV1
Abbreviation	Recombinant Cynomolgus monkey FASLG protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P63308
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	QLFHLQKELAELESTSQKHTASSLEKQIGHPSPPPEKKEQRKVAHLTGKPNRSMPLEWEDTYGIVLLSGVKYKKGGLVINETGLYFVYSKVYFRGQSCTNLPLSHKVYMRNSKYPQDLVMMEGKMSYCTTGQMWAHSSYLGAVFNLTSDHLYVNVSELSLVNFEESQTFFGLYKL
Research Area	Cell Biology
Source	E.coli
Target Names	FASLG
Protein Names	Recommended name: Tumor necrosis factor ligand superfamily member 6 Alternative name(s): CD95 ligand Short name= CD95-L Fas antigen ligand Short name= Fas ligand Short name= FasL CD_antigen= CD178 Cleaved into the f
Expression Region	102-280aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	26.5 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.