




Recombinant Human Acetylcholine receptor subunit alpha (CHRNA1), partial

Product Code	CSB-YP005386HUe9
Abbreviation	Recombinant Human CHRNA1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P02708
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SEHETRLVAKLFKDYSSVVRPVEDHRQVVEVTVGLQLIQLINVDEVNQIVTTNV RLKQGDMVDLPRPSCVTLGVPLFSHLQNEQWVDYNLKWNPDDYGGVKKIHIP SEKIWRPDLVLYNNADGDFAIVKFTKVLLQYTGHTWTPPAIFKSYCEIIVTHFPF DEQNCSMKLGTWTYDGSVVAINPESDQPDLSNFMESGEWVIKESRGWKHSV TYSCCPDTPYLDITYHFVMQRL
Research Area	Neuroscience
Source	Yeast
Target Names	CHRNA1
Protein Names	Recommended name: Acetylcholine receptor subunit alpha
Expression Region	21-255aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-PDI-tagged
Mol. Weight	86.0 kDa
Protein Length	Partial of Isoform 2
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.