




Recombinant Sandfly fever sicilian virus

Nucleoprotein (N)

Product Code	CSB-EP321019SCE
Abbreviation	Recombinant Sandfly fever sicilian virus N protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12793
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Sandfly fever sicilian virus (SFS)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MDEYQKIAVEFGGEQAIDETVIQDWLQAFAYQGFDARTIIHNLVQLGGKSWEED AKKMIILSLTRGNKPKKMVERMSPEGAREVKSLVAKYKIVEGRPGRNGITLSRV LQPWLGGQSKLWKWLKTSYQSQGAQWTALCGQTYPRQMMHPSFAGLIDPS LDQEDFNAVLD AHLFLFMFSKTINVSLRGAQKRDIEESFSQPMLAAINSSFIDN TQRR AFLTKFGILTSGARATAVVKKIAEVYRKLE
Research Area	Signal Transduction
Source	E.coli
Target Names	N
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= Protein N
Expression Region	1-246aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	32.0 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.