




Recombinant Mouse Aryl hydrocarbon receptor nuclear translocator-like protein 1 (Arntl), partial

Product Code	CSB-EP895734MO1f3
Abbreviation	Recombinant Mouse Arntl protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WTL8
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	EYAEHQGRIKNAREAHSQIEKRRRDKMNSFIDELASLVPTCNAMSRKLDKLTV LRMAVQHMKTLRGATNPYTEANYKPTFLSDDELKHLILRAADGFLFVVGCDRG KILFVSESVFKILNYSQNDLIGQSLFDYLHPKDIKVKEQLSSSDTAPRERLIDAK TGLPVKTDITPGPSRLCSGARRSFFCRMKCNRPVSVKVEDKDFASTCSKKKAD RKSFCTIHSTGYLKSWPPTKMGLDEDNEPDNEGCLNSCLVAIGRLHSHMVPQ PANGEIRVKSMEYVSRHAIDGKFVVDQRATAILAYLPQELLGTSCYEYFHQDD IGHLAECHRQVLQTREKITTNCYKFKIKDGSFITLRSRWFSFMNPWTKEVEYIV STNTVVLANVLE
Research Area	Cardiovascular
Source	E.coli
Target Names	Arntl
Expression Region	49-434aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-Myc-tagged
Mol. Weight	46.5 kDa
Protein Length	Partial of Isoform 3
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.