




# Recombinant Human Collagen alpha-6 (VI) chain (COL6A6), partial

<b>Product Code</b>	CSB-EP409861HUb0
<b>Abbreviation</b>	Recombinant Human COL6A6 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A6NMZ7
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	<p>VDTEEADIYLLIDGSGSTQATDFHEMKTFLSEVVGMFNIAPHKVRVGAVQYAD            SWDLEFEINKYSNKQDLGKAIENIRQMGGNTNTGAALNFTLSLLQKAKKQRGN            KVPCHLVVLTNGMSKDSILEPANRLREEHIRVYAIGIKEANQTQLREIAGEEKRV            YYVHDFDALKDIRNQVVQEICTEEACKEMKADIMFLVDSSGSIGPENFSKMKTF            MKNLVSKSQIGPDRVQIGVVQFSDINKEEFQLNRFMSQSDISNAIDQMAHIGQT            TLTGSALSFVSQYFSPTKGARNIRKFLILITDGEAQDIVKEPAVVLRRQEGVIIYS            VGVFGSNVTQLEEISGRPEMVFYVENFDILQRIEDDLVFGICSPREECKRIEVL            VVVIDSSGSIDYDEYNIMKDFMIGLVKKADVGKNQVRFGALKYADDPEVLVYL            DDFGKLEVISVLQNDQAMGGSTYAEALGFSDHMFTEARG SRLNKGVPQVLI            VITDGESHADKLNATAKALRDKGILVLA VGIDGANPVELLAMAGSSDKYFFVE            TFGGLKGFSDV</p>
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	COL6A6
<b>Expression Region</b>	430-982aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	67.4kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.