




Recombinant Escherichia coli Prophage outer membrane lipoprotein RzoR (rzoR)

Product Code	CSB-EP349921ENV
Abbreviation	Recombinant E.coli rzoR protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P58042
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli (strain K12)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	CTSKQSVSQCVCVKKPPPPAWIMQPPPDWQTPLNGIISPSGNDW
Research Area	Others
Source	E.coli
Target Names	rzoR
Protein Names	Recommended name: Outer membrane lipoprotein Rz1 from lambdoid prophage Rac
Expression Region	20-61aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	19.9 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.