




# Recombinant Vaccinia virus Protein B5 (PS/HR), partial

<b>Product Code</b>	CSB-EP321891VAA1
<b>Abbreviation</b>	Recombinant Vaccinia virus PS/HR protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P21115
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Vaccinia virus (strain Copenhagen) (VACV)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	YSTCTVPTMNNAKLTSTETSFNNNQKVTFTCDQGYHSSDPNAVCECDKWKYE NPCKKMCTVSDYISELYNKPLYEVNSTMTLSCNGETKYFRCEEKNGNTSWND TVTCPNAECQPLQLEHGSCQPVKEKYSFGYEYMTINCDVGYEVIGASYISCTAN SWNVIPSCQQKCDIPSLSNGLISGSTFSIGGVIHLCKSGFILTGSPSSTCIDGK WNPVLPICVRTNEEFDPVDDGPDDDETDLKLSKDVVQYEQEIESLEATYH
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	PS/HR
<b>Protein Names</b>	Recommended name: Plaque-size/host range protein Alternative name(s): Protein B5
<b>Expression Region</b>	18-279aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	36.5 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.